

Pe Rabbit Anti-Mouse Cd268/Baff-R Monoclonal Antibody For Flow Cytometry - Q9D8D0

Item Number: CM0008084

Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse CD268/BAFF-R (B-cell maturation defect/B-cell-activating factor receptor). Validated for flow cytometry, suitable for B-cell receptor signaling and mature B-cell survival studies.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD268/BAFF-R mAb
Remarks/Alias	Bcmd; Baffr; Bcmd1; BAFF-R; Bcmd-1; Lvis22; 2010006P15Rik
Species	Mouse
GeneID (Human)	115650
GeneID	72049
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 10-71 of mouse CD268/BAFF-R (NP_082351.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9D8D0
Protein Weight	19kDa
Shipping	Ice bag