

Pe Rabbit Anti-Mouse Cd272/Btla Monoclonal Antibody For Fc - Q7Tsa3

Item Number: CM0009654



Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse CD272/BTLA, validated for flow cytometry. B- and T-lymphocyte attenuator involved in immune regulation. Ideal for immunological research applications.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD272/BTLA mAb
Remarks/Alias	A630002H24
Species	Mouse
GeneID (Human)	151888
GeneID	208154
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 30-176 of mouse CD272/BTLA (NP_001032808.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q7TSA3
Protein Weight	22kDa/34kDa
Shipping	Ice bag