


Pe Rabbit Anti-Mouse Cd274/Pd-L1 Monoclonal Antibody For Flow Cytometry - Q9Ep73

Item Number: CM0011843



Introduction

PE-conjugated rabbit monoclonal antibody against mouse PD-L1 (CD274), targeting the extracellular domain. Validated for flow cytometry in mouse samples; suitable for immune checkpoint and tumor immunology research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD274/PD-L1 mAb
Remarks/Alias	B7h1, Pdl1, Pdcd1l1, Pdcd1lg1, A530045L16Rik
Species	Mouse
GeneID (Human)	29126
GeneID	60533
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9EP73
Protein Weight	32kDa
Shipping	Ice bag