

Pe Rabbit Anti-Mouse Cd28 Monoclonal Antibody For Flow Cytometry - P31041

Item Number: CM0004635



Introduction

PE-conjugated rabbit monoclonal antibody against mouse CD28, validated for flow cytometry. Ideal for studying T-cell activation, costimulation, and immune homeostasis. Reacts with mouse.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD28 mAb
Remarks/Alias	Cd28
Species	Mouse
Gene ID (Human)	940
Gene ID	12487
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 20-149 of mouse CD28 (NP_031668.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P31041
Protein Weight	25kDa
Shipping	Ice bag