

Pe Rabbit Anti-Mouse Cd301 (Mgl1/2) Mab For Flow Cytometry - P49300 / Q8Jzn1

Item Number: CM0012697

Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse CD301 (MGL1/2) for flow cytometry. Recognizes CLEC10A lectin on macrophages; aids research in tumor immunity, adipose inflammation, and infection.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD301 (MGL1/2) mAb
Remarks/Alias	Mgl; Mgl1; CD301a; M-ASGP-BP-1
Species	Mouse
GenelD (Human)	10462/54551
GenelD (Mouse)	17312/216864
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 57-304 of mouse CD301a (NP_034926.1))
Source	Rabbit
Category	Monoclonal Antibodies
Application	Flow Cytometry
Cross Reactivity	Mouse
SWISS	P49300/Q8JZN1
Shipping	Ice bag