

Pe Rabbit Anti-Mouse Cd301B/Mgl2 Monoclonal Antibody For Fc - Q8Jzn1

Item Number: CM0008665

Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse CD301b (MGL2), a C-type lectin membrane glycoprotein. Validated for flow cytometry (FC) applications. Ideal for immunology and glycobiology research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD301b/MGL2 mAb
Remarks/Alias	CD301b
Species	Mouse
GeneID (Human)	10462
GeneID	216864
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 233-332 of mouse CD301b/MGL2 (NP_660119.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q8JZN1
Protein Weight	38 kDa
Shipping	Ice bag