


Pe Rabbit Anti-Mouse Cd304/Neuropilin-1 Monoclonal Antibody For Fc - P97333

Item Number: CM0011390



Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse CD304/Neuropilin-1, validated for flow cytometry (FC) applications. Suitable for investigating neurogenesis, angiogenesis, and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD304/Neuropilin-1 mAb
Remarks/Alias	Nrp; NP-1; Npn1; NPN-1; C530029I03
Species	Mouse
GeneID (Human)	8829
GeneID	18186
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-856 of mouse CD304/Neuropilin-1 (NP_032763.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P97333
Protein Weight	103kDa
Shipping	Ice bag