

Pe Rabbit Anti-Mouse Cd32 Monoclonal Antibody For Fc - P08101

Item Number: CM0012583

Introduction

PE-conjugated rabbit monoclonal antibody against mouse CD32 (FcγRIIB), validated for flow cytometry. Suitable for studying IgG receptor function in lymphoid and myeloid cells. Uniprot: P08101.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD32 mAb
Remarks/Alias	CD32, Fcgr2, Fcr-2, Fcr-3, Ly-17, LyM-1, Lym-1, fcRII, FcgRII, Fcgr2a, Ly-m20, Fc[g]RII, F630109E10Rik
Species	Mouse
Gene ID (Human)	2213
Gene ID	14130
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-207 of mouse CD32 (NP_034317.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P08101
Protein Weight	31kDa/33kDa/36kDa
Shipping	Ice bag