

Pe Rabbit Anti-Mouse Cd34 Monoclonal Antibody For Fc - Q64314

Item Number: CM0010280

Introduction

PE-labeled rabbit monoclonal antibody against mouse CD34, validated for flow cytometry (FC). Suitable for studying hematopoietic stem/progenitor cell adhesion and early hematopoiesis. Reacts with mouse CD34. UniProt: Q64314.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD34 mAb
Remarks/Alias	CD34;CD34 molecule;GIG3;MORT1
Species	Mouse
GeneID (Human)	947
GeneID	12490
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-287 of mouse CD34 (NP_598415.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q64314
Protein Weight	41kDa
Shipping	Ice bag