

Pe Rabbit Anti-Mouse Cd35 Mab - Cd35

Item Number: CM0009516

Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse CD35 (Cr1, C3DR) for flow cytometry (FC). Reacts with mouse, protein weight 113 kDa. Ideal for immunology research on complement receptors.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD35 mAb
Remarks/Alias	Cr1; C3DR; CD21; CD35; Cr-1; Cr-2
Species	Mouse
Gene ID (Human)	1380
Gene ID	12902
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 21-396 of mouse CD35 (NP_001355694.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
Protein Weight	113 kDa
Shipping	Ice bag