

Pe Rabbit Anti-Mouse Cd352/Slamf6 Monoclonal Antibody For Flow Cytometry - Q9Et39

Item Number: CM0008512

Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse CD352/SLAMF6 (Ly108/NTB-A) validated for flow cytometry. Suitable for studying SLAM family receptor-mediated immune cell signaling and modulation in mouse hematopoietic samples.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD352/SLAMF6 mAb
Remarks/Alias	KAL1; NTBA; KAL1b; Ly108; NTB-A; SF2000
Species	Mouse
GeneID (Human)	114836
GeneID	30925
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-235 of mouse CD352/SLAMF6 (NP_001334115.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9ET39
Protein Weight	39kDa
Shipping	Ice bag