

Pe Rabbit Anti-Mouse Cd366 (Tim-3) Monoclonal Antibody For Flow Cytometry - Q8Vim0

Item Number: CM0012169



Introduction

PE-conjugated rabbit monoclonal antibody against mouse CD366 (TIM-3), validated for flow cytometry. Specific for mouse, expressed on Th1, CD8+ T cells, macrophages, and dendritic cells. Ideal for immune checkpoint research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD366 (TIM-3) mAb
Remarks/Alias	Tim3; TIM-3; Timd3
Species	Mouse
GeneID (Human)	84868
GeneID	171285
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-191 of mouse CD366 (TIM-3) (NP_599011.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q8VIM0
Protein Weight	22kDa/31kDa
Shipping	Ice bag