

Pe Rabbit Anti-Mouse Cd40 Monoclonal Antibody For Flow Cytometry - P27512

Item Number: CM0009238

Introduction

PE-conjugated rabbit monoclonal antibody specific to mouse CD40 (TNFRSF5, Bp50), validated for flow cytometry. Targets extracellular domain (aa 20-189). Ideal for B cell activation and immune signaling studies.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD40 mAb
Aliases	IGM; p50; Bp50; GP39; IMD3; TRAP; HIGM1; T-BAM; Tnfrsf5
Species	Mouse
Gene ID (Human)	958
Gene ID	21939
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 20-189 of mouse CD40 (NP_035741.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P27512
Protein Weight	32kDa
Shipping	Ice bag