

Pe Rabbit Anti-Mouse Cd43 Mab - P15702

Item Number: CM0005509



Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse CD43 (Leukosialin). Validated for flow cytometry, suitable for immunophenotyping, T-cell activation, and migration studies in mouse models.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD43 mAb
Remarks/Alias	Cd43; Ly48; Galgp; Ly-48; A630014B01Rik
Species	Mouse
GeneID (Human)	6693
GeneID	20737
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-248 of mouse CD43 (NP_033285.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P15702
Protein Weight	40 kDa
Shipping	Ice bag