


Pe Rabbit Anti-Mouse Cd44 Monoclonal Antibody For Fc - P15379

Item Number: CM0009479



Introduction

PE-conjugated rabbit monoclonal antibody against mouse CD44 (Ly-24, Pgp-1, HERMES). Validated for flow cytometry (FC). Suitable for studying cell adhesion, lymphocyte homing, and hyaluronan binding.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD44 mAb
Remarks/Alias	Ly-24; Pgp-1; HERMES
Species	Mouse
Gene ID (Human)	960
Gene ID	12505
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 23-176 of mouse CD44 (NP_033981.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P15379
Protein Weight	85kDa
Shipping	Ice bag