

Pe Rabbit Anti-Mouse Cd45 Mab For Fc - P06800

Item Number: CM0009683

Introduction

Reliable PE-conjugated rabbit monoclonal antibody specific to mouse CD45 (PTPRC, Leukocyte common antigen). Optimized for flow cytometry (FC), ideal for immunophenotyping and diagnostic raw material sourcing.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD45 mAb
Remarks/Alias	loc; B220; Cd45; L-CA; Ly-5; T200; CD45R; Lyt-4
Species	Mouse
GeneID (Human)	5788
GeneID	19264
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-566 of mouse CD45 (NP_001104786.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P06800
Protein Weight	130kDa/140kDa/144kDa
Shipping	Ice bag