

Pe Rabbit Anti-Mouse Cd47 Monoclonal Antibody For Fc - Q61735

Item Number: CM0010063

Introduction

PE-conjugated rabbit monoclonal antibody against mouse CD47 (IAP) for flow cytometry. Validated for mouse samples, targeting the leukocyte surface antigen involved in cell adhesion and immune regulation. UniProt Q61735.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD47 mAb
Remarks/Alias	IAP; Itgp; 9130415E20Rik; B430305P08Rik
Species	Mouse
GeneID (human)	961
GeneID	16423
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-135 of mouse CD47 (NP_034711.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q61735
Protein Weight	35kDa
Shipping	Ice bag