

Pe Rabbit Anti-Mouse Cd5 Monoclonal Antibody For Fc - P13379

Item Number: CM0011559

Introduction

PE-conjugated rabbit monoclonal antibody specific for mouse CD5 (Ly-1 antigen), validated for flow cytometry. Suitable for T-cell and B1a cell research, immune tolerance, and signaling studies.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD5 mAb
Remarks/Alias	Ly-1; Ly-A; Ly-12; Lyt-1
Species	Mouse
GeneID (Human)	921
GeneID	12507
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 25-371 of mouse CD5 (NP_031676.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P13379
Protein Weight	54kDa
Shipping	Ice bag