

Pe Rabbit Anti-Mouse Cd50/Icam-3 (Icam-5/Telencephalin) Monoclonal Antibody For Fc - Q60625

Item Number: CM0006018

Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse CD50/ICAM-3 (ICAM-5/Telencephalin), validated for flow cytometry. Ideal for neuronal adhesion and neuroimmunology research from CamelBio.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD50/ICAM-3 mAb
Remarks/Alias	TLN; CD50; Tlcn; Icam3; ICAM-3; ICAM-5
Species	Mouse
GeneID (Human)	7087
GeneID (Mouse)	15898
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-227 of mouse CD50/ICAM-3 (NP_032345.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q60625
Protein Weight	97kDa
Shipping	Ice bag