

Pe Rabbit Anti-Mouse Cd6 Mab - Q61003

Item Number: CM0004942

Introduction

PE Rabbit anti-Mouse CD6 monoclonal antibody for flow cytometry. Detects mouse CD6, a T-cell antigen involved in adhesion, costimulation, and innate immunity. MW: 64/68/72 kDa. Swiss-Prot: Q61003. Ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD6 mAb
Remarks/Alias	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Species	Mouse
GeneID (Human)	923
GeneID	12511
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-155 of mouse CD6 (NP_033982.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q61003
Protein Weight	64kDa/68kDa/72kDa
Shipping	Ice bag