

Pe Rabbit Anti-Mouse Cd62E Monoclonal Antibody For Fc - Q00690

Item Number: CM0010851

Introduction

PE-conjugated rabbit monoclonal antibody against mouse CD62E (E-selectin), validated for flow cytometry. Suitable for immunoadhesion and leukocyte trafficking studies in murine models.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD62E mAb
Remarks/Alias	Elam; CD62E; ELAM-1; LECAM2; E-selectin
Species	Mouse
GeneID (human)	6401
GeneID	20339
Immunogen	Recombinant protein corresponding to a sequence within amino acids 22-557 of mouse CD62E (NP_035475.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q00690
Protein Weight	66kDa
Shipping	Ice bag