

Pe Rabbit Anti-Mouse Cd62L Monoclonal Antibody For Flow Cytometry - P18337

Item Number: CM0009601

Introduction

PE Rabbit anti-Mouse CD62L monoclonal antibody for flow cytometry. Recognizes L-selectin (CD62L), validated in mouse. Suitable for studying lymphocyte homing and leukocyte rolling. SWISS: P18337.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD62L mAb
Remarks/Alias	Lnhr; CD62L; LAM-1; Ly-22; Lyam1; LECAM1; Ly-m22; Lyam-1; LECAM-1; L-selectin
Species	Mouse
GenelD (human)	6402
GenelD	20343
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 39-332 of mouse CD62L (NP_035476.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P18337
Protein Weight	42kDa
Shipping	Ice bag