

Pe Rabbit Anti-Mouse Cd68 Monoclonal Antibody For Flow Cytometry (Intracellular) - P31996

Item Number: CM0012762

Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse CD68 (macrosialin) for intracellular flow cytometry applications. Suitable for mouse macrophage and dendritic cell staining. High specificity and bright PE signal.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD68 mAb
Remarks/Alias	Lamp4; gp110; Scard1
Species	Mouse
GeneID (Human)	968
GeneID	12514
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-291 of mouse CD68 (NP_001277987.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	P31996
Protein Weight	35kDa
Shipping	Ice bag