

Pe Rabbit Anti-Mouse Cd70 Mab For Fc - O55237

Item Number: CM0008570

Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse CD70. Validated for flow cytometry (FC). Used for investigating CD70-CD27 signaling in T cell activation and B cell immune surveillance.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD70 mAb
Remarks/Alias	Cd27l; CD27LG; Tnfsf7; Tnlg8a
Species	Mouse
GeneID (Human)	970
GeneID	21948
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 47-195 of mouse CD70 (NP_035747.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	O55237
Protein Weight	22 kDa
Shipping	Ice bag