

## Pe Rabbit Anti-Mouse Cd79B Monoclonal Antibody For Fc - P15530

Item Number: CM0010048



### Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse CD79b (B29, Ig-beta, B-cell-specific glycoprotein B29), validated for flow cytometry (FC) applications. Ideal for studying B-cell receptor signaling and adaptive immunity.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | PE Rabbit anti-Mouse CD79b mAb                                                                                    |
| Remarks/Alias    | B29; Igb; Igbeta; Ig-beta                                                                                         |
| Species          | Mouse                                                                                                             |
| GeneID (Human)   | 974                                                                                                               |
| GeneID           | 15985                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 26-158 of mouse CD79b (NP_032365.1) |
| Source           | Rabbit                                                                                                            |
| Category         | Monoclonal Antibodies                                                                                             |
| Application      | FC                                                                                                                |
| Cross Reactivity | Mouse                                                                                                             |
| SWISS            | P15530                                                                                                            |
| Protein Weight   | 25kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |