

Pe Rabbit Anti-Mouse Cd80 Mab For Fc - Q00609

Item Number: CM0009522



Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse CD80, a key T-cell costimulatory molecule. Validated for flow cytometry, ideal for immunology research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD80 mAb
Remarks/Alias	B71; Ly53; TSA1; Cd28I; Ly-53; MIC17
Species	Mouse
GeneID (Human)	941
GeneID	12519
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 38-245 of mouse CD80, NP_033985.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q00609
Protein Weight	34kDa
Shipping	Ice bag