

Pe Rabbit Anti-Mouse Cd82 Monoclonal Antibody For Flow Cytometry - P40237

Item Number: CM0004410

Introduction

PE-conjugated rabbit anti-mouse CD82 monoclonal antibody for flow cytometry. Recognizes mouse CD82, a tetraspanin involved in immune regulation, adhesion, and angiogenesis. Suitable for mouse splenocyte and kidney cell analysis.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD82 mAb
Remarks/Alias	C33; IA4; Kai1; Tspan27
Species	Mouse
GeneID (human)	3732
GeneID	12521
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-266 of mouse CD82 (NP_001258390.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P40237
Protein Weight	30kDa
Shipping	Ice bag