

Pe Rabbit Anti-Mouse Cd83 Monoclonal Antibody For Flow Cytometry - O88324

Item Number: CM0010173

Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse CD83, validated for flow cytometry (FC). Ideal for immune cell research, dendritic cell and lymphocyte studies, and immune response modulation investigations.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD83 mAb
Remarks/Alias	mCD83
Species	Mouse
GenelD (Human)	9308
GenelD	12522
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 22-133 of mouse CD83 (NP_033986.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	O88324
Protein Weight	21KDa
Shipping	Ice bag