

Pe Rabbit Anti-Mouse Cd84 Mab - Q18Pi6

Item Number: CM0009651

Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse CD84 (SLAMF5) for flow cytometry. Ideal for studies on SLAM family receptors, adaptive immunity, and NK/T cell function.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD84 mAb
Remarks/Alias	CDw84; SLAMF5; A130013D22Rik
Species	Mouse
GeneID (Human)	8832
GeneID	12523
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-221 of mouse CD84 (NP_038517.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q18Pi6
Protein Weight	37kDa
Shipping	Ice bag