

Pe Rabbit Anti-Mouse Cd86 Mab - P42082

Item Number: CM0010055

Introduction

PE-conjugated rabbit monoclonal antibody against mouse CD86, validated for flow cytometry. Suitable for T-cell activation and costimulation pathway studies.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD86 mAb
Remarks/Alias	B7; B70; MB7; B7-2; B7.2; CLS1; Ly58; ETC-1; Ly-58; MB7-2; Cd28l2; TS/A-2; CD86
Species	Mouse
Gene ID (Human)	942
Gene ID	12524
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P42082
Protein Weight	35 kDa
Shipping	Ice bag