

Pe Rabbit Anti-Mouse Cd88/C5Ar Monoclonal Antibody For Flow Cytometry - P30993

Item Number: CM0009528



Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse CD88/C5aR, validated for flow cytometry (FC). Reacts with mouse, predicted protein weight 39 kDa. Suitable for immunology and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD88/C5aR mAb
Remarks/Alias	C5aR; C5r1; Cd88; D7Msu1
Species	Mouse
GeneID (Human)	728
GeneID (Mouse)	12273
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-37 of mouse CD88/C5aR (NP_031603.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P30993
Protein Weight	39kDa
Shipping	Ice bag