

Pe Rabbit Anti-Mouse Cd90.2 Monoclonal Antibody For Flow Cytometry - P01831

Item Number: CM0010641

Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse CD90.2 (Thy-1) for flow cytometry. Suitable for mouse samples; aids in immunophenotyping and cell-cell interaction studies.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD90.2 mAb
Remarks/Alias	T25; CD90; Thy-1; Thy1.1; Thy1.2; Thy-1.2
Species	Mouse
Gene ID (Human)	7070
Gene ID	21838
Immunogen	Recombinant protein of mouse CD90.2
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01831
Protein Weight	18kDa
Shipping	Ice bag