


Pe Rabbit Anti-Mouse Cd95/Fas Mab - P25446

Item Number: CM0012195



Introduction

PE-conjugated rabbit monoclonal antibody against mouse CD95/FAS for flow cytometry. Recognizes mouse tumor necrosis factor receptor superfamily member 6 (SWISS P25446); suitable for apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD95/FAS mAb
Remarks/Alias	lpr; APO1; APT1; CD95; TNFR6; Tnfrsf6
Species	Mouse
Gene ID (Human)	355
Gene ID	14102
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-169 of mouse CD95/FAS (NP_032013.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P25446
Protein Weight	37 KDa
Shipping	Ice bag