


Pe Rabbit Anti-Mouse Cd96/Tactile Monoclonal Antibody For Flow Cytometry - Q3U0X8

Item Number: CM0012427



Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse CD96/TACTILE, validated for flow cytometry. Ideal for studying T cell and NK cell adhesion and late-phase immune responses. Q3U0X8.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD96/TACTILE mAb
Remarks/Alias	Tactile; 1700109I12Rik
Species	Mouse
Gene ID (Human)	10225
Gene ID	84544
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-536 of mouse CD96/TACTILE (NP_115854.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q3U0X8
Protein Weight	67kDa
Shipping	Ice bag