


Pe Rabbit Anti-Mouse E-Cadherin/Cd324 Monoclonal Antibody For Fc - P09803

Item Number: CM0008566



Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse E-Cadherin/CD324 (SWISS: P09803). Validated for flow cytometry (FC). Suitable for studying cell adhesion, epithelial biology, and cancer research.

[Learn More](#)

Parameter	Specification
Product Name	PE Rabbit anti-Mouse E-Cadherin/CD324 mAb
Remarks/Alias	Um; UVO; Ecad; ARC-1; E-cad; L-CAM
Species	Mouse
GeneID (Human)	999
GeneID	12550
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 157-369 of mouse E-Cadherin/CD324 (NP_033994.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P09803
Protein Weight	98kDa
Shipping	Ice bag