

Pe Rabbit Anti-Mouse Foxp3 Monoclonal Antibody For Fc (Intra) - Q99Jb6

Item Number: CM0004119



Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse FOXP3 (scurfin), a key transcriptional regulator in regulatory T-cells (Tregs). Validated for intracellular flow cytometry, this antibody enables detection and analysis of Treg populations in mouse samples.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse FOXP3 mAb
Remarks/Alias	sf; JM2; scurfin
Species	Mouse
Gene ID (Human)	50943
Gene ID	20371
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of mouse Foxp3 (NP_473380.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	Q99JB6
Protein Weight	47kDa
Shipping	Ice bag