


Pe Rabbit Anti-Mouse I-A Monoclonal Antibody For Fc - P14483

Item Number: CM0010672



Introduction

PE-conjugated rabbit monoclonal antibody specific for mouse I-A MHC class II alloantigens (H-2I, I-A, I-E). Validated for flow cytometry. Suitable for immune cell phenotyping and immunological research.

[Learn More](#)

Parameter	Specification
Product Name	PE Rabbit anti-Mouse I-A mAb
Remarks/Alias	H-2I; I-Ad, -Ed; I-Eb and I-Ek MHC class II alloantigens; Ia Ag; MHC II
Species	Mouse
GeneID (human)	474382
GeneID	14961
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P14483
Protein Weight	30kDa
Shipping	Ice bag