

Pe Rabbit Anti-Mouse Icos/Cd278 Monoclonal Antibody For Fc - Q9Wvs0

Item Number: CM0011122

Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse ICOS/CD278, validated for flow cytometry (FC). Recognizes the inducible T-cell costimulator, critical for T-cell activation and T-B cell interactions, in mouse samples.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse ICOS/CD278 mAb
Remarks/Alias	H4; CCLP; AILIM; CRP-1; Ly115
Species	Mouse
GeneID (Human)	29851
GeneID	54167
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-142 of Mouse ICOS/CD278 (NP_059508.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9WVS0
Protein Weight	22kDa
Shipping	Ice bag