

Pe Rabbit Anti-Mouse Iga (Fc) Monoclonal Antibody

Item Number: CM0012556

Introduction

PE-conjugated rabbit monoclonal antibody specific for mouse IgA Fc region. Designed for intracellular flow cytometry in mouse samples. Sourced from rabbit. Shipped on ice. Ideal for immunological research and assay development.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse IgA (Fc) mAb
Species	Mouse
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 99-344 of mouse IgA
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (Intracellular)
Cross Reactivity	Mouse
Shipping	Ice bag