


Pe Rabbit Anti-Mouse IgG1 Monoclonal Antibody For Fc - P01868

Item Number: CM0005101



Introduction

PE-conjugated rabbit monoclonal antibody against mouse IgG1, suitable for flow cytometry (FC). This antibody is specific for the IgG1 isotype and can be used for mouse immunoglobulin profiling and B cell analysis.

[Learn More](#)

Parameter	Description
Product Name	PE Rabbit anti-Mouse IgG1 mAb
Remarks/Alias	IgG1
Species	Mouse
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 98-324 of mouse IgG1 (P01868)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01868
Protein Weight	36kDa
Shipping	Ice bag