

Pe Rabbit Anti-Mouse Igg2B Mab For Fc - P01867-2

Item Number: CM0008323

Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse IgG2b constant region (isoform 2). Validated for flow cytometry with cross-reactivity to mouse. Suitable for immunoglobulin isotyping and immune profiling studies.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse IgG2b mAb
Remarks/Alias	IgG2b; Igh-3; gamma2b
Species	Mouse
GeneID	16016
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 14-113 of mouse IgG2b (P01867-2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01867-2
Protein Weight	37kDa
Shipping	Ice bag