

Pe Rabbit Anti-Mouse IgM Monoclonal Antibody For Fc - P01872

Item Number: CM0008613



Introduction

PE-conjugated rabbit monoclonal antibody against mouse IgM constant region (P01872). Validated for flow cytometry. Reacts specifically with mouse IgM, enabling B cell detection and humoral immunity studies.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse IgM mAb
Remarks/Alias	Igm; muH; Igh6; Igh-6; Igh-M
Species	Mouse
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 105-454 of Mouse IgM (P01872)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01872
Protein Weight	50kDa/53kDa
Shipping	Ice bag