

Pe Rabbit Anti-Mouse IL-13 Monoclonal Antibody For Flow Cytometry (Intracellular) - P20109

Item Number: CM0011563

Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse IL-13, suitable for intracellular flow cytometry. Ideal for studies in allergic inflammation, immune regulation, and cytokine signaling pathways. - P20109

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse IL-13 mAb
Remarks/Alias	IL-13
Species	Mouse
GenelD (Human)	3596
GenelD	16163
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-131 of mouse IL-13 (NP_032381.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	P20109
Protein Weight	14kDa
Shipping	Ice bag