

Pe Rabbit Anti-Mouse IL-17F Monoclonal Antibody - Q7Tni7

Item Number: CM0006073

Introduction

PE-conjugated rabbit monoclonal antibody against mouse IL-17F, designed for intracellular flow cytometry. Explore Th17 cell biology, antimicrobial immunity, and inflammatory responses with this validated reagent.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse IL-17F mAb
Remarks/Alias	IL-17F
Species	Mouse
GeneID (Human)	112744
GeneID	257630
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-161 of mouse IL-17F (NP_665855.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	Q7TNI7
Protein Weight	18 kDa
Shipping	Ice bag