


Pe Rabbit Anti-Mouse IL-1B Monoclonal Antibody For Flow Cytometry (Intracellular) - P10749

Item Number: CM0005190



Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse IL-1 β (P10749) for intracellular flow cytometry. Validated for mouse samples. Ideal for immune cell activation, inflammation, and pyroptosis research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse IL-1 β mAb
Remarks/Alias	IL-1b; IL-1beta
Species	Mouse
GeneID (human)	3553
GeneID	16176
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 118-269 of mouse IL1 β (NP_032387.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	P10749
Protein Weight	31 kDa
Shipping	Ice bag