

Pe Rabbit Anti-Mouse IL-22 Monoclonal Antibody For Intracellular Flow Cytometry - Q9Jjy9

Item Number: CM0010723



Introduction

PE-conjugated rabbit monoclonal antibody against mouse IL-22 for intracellular flow cytometry. Reacts specifically with mouse IL-22. Ideal for studying IL-22-mediated tissue regeneration, barrier function, and inflammatory responses.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse IL22 mAb
Remarks/Alias	IL-22; If2b1; Iltif; IL-22a; ILTIFa
Species	Mouse
GeneID (human)	50616
GeneID	50929
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	Q9JJY9
Protein Weight	20kDa
Shipping	Ice bag