


Pe Rabbit Anti-Mouse IL-33R/St2 Mab - P14719

Item Number: CM0005883



Introduction

PE-conjugated rabbit monoclonal antibody against mouse IL-33R/ST2 (ST2). Suitable for flow cytometry (FC) in immunological studies of IL-33 signaling and macrophage polarization. Cross-reacts with mouse.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse IL-33R/ST2 mAb
Remarks/Alias	T1; St2; DER4; Ly84; ST2L; Fit-1; T1/ST2; St2-rs1
Species	Mouse
Gene ID (Human)	9173
Gene ID	17082
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-332 of mouse IL-33R/ST2 (NP_001020773.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P14719
Protein Weight	38kDa/65kDa
Shipping	Ice bag