

Pe Rabbit Anti-Mouse Il-9 Monoclonal Antibody For Flow Cytometry (Intracellular) - P15247

Item Number: CM0005775



Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse IL-9 (P40). Suitable for intracellular flow cytometry. Reacts specifically with mouse IL-9, a key cytokine in Th2 responses and mast cell biology. Ideal for immune function studies.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse IL-9 mAb
Remarks/Alias	P40; IL-9
Species	Mouse
GeneID (Human)	3578
GeneID	16198
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-114 of mouse IL-9 (NP_032399.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	P15247
Protein Weight	16kDa
Shipping	Ice bag