

## Pe Rabbit Anti-Mouse IL1R1/Cd121A Monoclonal Antibody For Flow Cytometry - P13504

Item Number: CM0012701

### Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse IL1R1/CD121a (P13504), validated for flow cytometry. Ideal for studying interleukin-1 receptor type 1 in inflammatory response pathways.

[Learn More](#)

| Parameter        | Value                                                                                                                    |
|------------------|--------------------------------------------------------------------------------------------------------------------------|
| Product Name     | PE Rabbit anti-Mouse IL1R1/CD121a mAb                                                                                    |
| Remarks/Alias    | IL-iR; CD121a; CD121b; IL-1R1; Il1r-1; IL-1R-1; IL-1RT1; IL-1RT-1; IL-1R-alpha                                           |
| Species          | Mouse                                                                                                                    |
| GenelD (Human)   | 3554                                                                                                                     |
| GenelD           | 16177                                                                                                                    |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 20-338 of mouse IL1R1/CD121a (NP_032388.1) |
| Source           | Rabbit                                                                                                                   |
| Category         | Monoclonal Antibodies                                                                                                    |
| Application      | FC                                                                                                                       |
| Cross Reactivity | Mouse                                                                                                                    |
| SWISS            | P13504                                                                                                                   |
| Protein Weight   | 67kDa                                                                                                                    |
| Shipping         | Ice bag                                                                                                                  |