

Pe Rabbit Anti-Mouse Il4 Monoclonal Antibody For Intracellular Fc - P07750

Item Number: CM0011876

Introduction

PE-labeled rabbit monoclonal antibody against mouse IL4 (BSF-1), validated for intracellular flow cytometry. Suitable for studying T-cell responses, inflammation, and neuroimmunology in mouse models.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse IL4 mAb
Remarks/Alias	IL-4; BSF-1
Species	Mouse
GenelD (Human)	3565
GenelD (Mouse)	16189
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-140 of mouse IL4 (NP_067258.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	P07750
Protein Weight	16 kDa
Shipping	Ice bag